




# CentriVap® Centrifugal Concentrators and Cold Traps



# Overview

- What is a Concentrator?
  - Evaporator that uses centrifuge motion, vacuum and optional heat to concentrate a liquid sample to dryness or a wet pellet.
- Applications that relate to current events
  - Paternity testing – Anna Nicole Smith trial
  - Drugs of abuse testing – Steroid use in athletes
  - Drug development – New drug discovery, diabetes, cancer, etc.






# Definition: Centrifugal Concentrator

A device that uses centrifugal force and vacuum (and optionally heat), to concentrate a selected solute from solution. Also called a vacuum concentrator.



## Definition: Centrifugal Force

- The component of apparent force on a body in curvilinear motion, as observed from that body, that is directed away from the center of curvature or axis of rotation.




## Definition: Low Vacuum Pressures

- Low vacuum pressures help vaporize solvents. The rate of evaporation of the solvent increases because there is a decrease in vapor pressure, which encourages the solvent molecules to pass from the liquid phase to the gaseous phase.




# Evaporation and Sublimation



- Evaporation occurs when a liquid transforms into a vapor
- Occasionally, at low pressure and temperature, a sample will freeze and undergo sublimation
- Sublimation occurs when a solid transforms directly to a vapor without going through the liquid phase


# Evaporation

- The rate of evaporation can be increased by reducing the system pressure and by heating the solvent.
- Low system pressure decreases the vapor pressure, and that encourages the solvent to go from liquid to gaseous phase.
- Increasing the temperature will increase the movement of the molecules, which also encourages the solvent molecules to enter the gaseous phase.
- By increasing the surface area, the rate of evaporation increases.



# Theory of Operation


- As the rotor turns, centrifugal force moves the solvent and solute to the bottom of the tube. This prevents foaming and/or bumping during the concentration process






# Theory of Operation

- By connecting a vacuum pump to the centrifugation chamber, the vapor pressure is reduced, increasing the speed of concentration.
- If the solvent is liquid, then it evaporates; if it is solid, sublimation occurs.




# Theory of Operation



Placing a cold trap in the vacuum line allows the user to collect the solvent.

# Theory of Operation

- A secondary chemical trap placed after the cold trap, physically traps acids, moisture, radioisotopes, or solvents to protect the vacuum pump.



# Solute Recovery

- The solvent is removed from the solute by evaporation and/or sublimation. Solvent recovery is 99%. The remaining sample pellet is nearly pure solute.



# History

- First manufactured in 1991
- Most recent redesign 2005
- New Products Refrigerated CentriVap and Acid Resistant CentriVap introduced 2006





# Overview of Products

# Features Common to All Concentrators







# CentriVap® Control Panel

Program - shows program from 1-9, R for run or S for stop.

Quick-Start™ Buttons initiate user-set Program 1, 2 or 3 with single push of button.

Heater Time - shows time selected for heat from 1-999 minutes.

Increase/Decrease Buttons - increase or decrease last selected set point parameter.



Temperature - shows heat (and Heat Boost, if activated) from OFF to 99° C or HI (100° C).

Preheat Button - turns on heater to preheat chamber prior to loading samples.

Run Time - shows time from 1-999 minutes. During run, time remaining also displays. Selecting "ON" allows CentriVap to run continuously.





# CentriVap® Centrifugal Concentrator

## Control Panel



### *Lightbulb feature!*

- Heat and run times are each set separately from 1 to 999 minutes so heat can be programmed to shut down, protecting samples from excessive heat exposure as they dry
- LCD display prompts user to set program parameters and displays current program



# CentriVap® Centrifugal Concentrator

## Control Panel



- Quick-Start™ One Button Start Up activates the rotor, heater, timers and vacuum pump with one button\*
- Preheat button activates the heater to begin elevating to set point temp
- Memory stores 9 user-set programs\*
- Audible alarm signals completion of set point run time



# CentriVap® Centrifugal Concentrator

## Common Features

- Lid latch with safety sensor
- Built-in vacuum delay prevents bumping
- Variety of rotors available
- Automatic vacuum release valve
- Epoxy-coated aluminum chamber, Teflon coated as an option





## CentriVap® Centrifugal Concentrator Common Features

- Heat Boost™ for extra heat controllable in 1° C increments through sidewalls of chamber\*
- Larger chamber
- Quick-Stop System
- Rear vacuum port



# CentriVap® Centrifugal Concentrator

***Lightbulb feature!***

**CentriZap™**



- Optional accessory allows viewing samples while rotor is turning. Samples appear to be standing still



- Plugs into outlet on back of CentriVap and holster attaches with two fasteners to side of cabinet

\*Exclusive feature

# CentriVap® Centrifugal Concentrator



## Accessories

- Glass lid provides additional protection from solvents and corrosives that attack acrylic
- Secondary Chemical Traps remove vapors and liquid that may be exhausted from the vacuum pump and into the laboratory
- Digital Electronic Vacuum Gauge





# CentriVap® DNA Centrifugal Concentrator



Incorporates a built-in diaphragm vacuum pump

- DNA rotor included for 0.5 – 2ml tubes
- Two 177 ml glass traps protect vacuum pump and environment\*
- Rotor accommodates up to 60 each 1.5 ml and 72 each 0.5 ml microcentrifuge tubes
- Accessory Microtiter Plate Rotor holds four standard 96-well plates or two deep well plates

***Lightbulb feature!***

\*Exclusive feature

**LABCONCO**



# CentriVap® DNA Centrifugal Concentrator

177 ml glass  
traps

Integrated  
diaphragm pump


DNA Rotor included

Control Panel

Larger chamber and  
rotors hold more  
samples for greater  
throughput

LCD Display

Microprocessor-  
controlled 300-watt  
heater



**LABCONCO**





# CentriVap® Centrifugal Concentration Systems


**Designed to rapidly concentrate  
multiple biological samples**



**LABCONCO®**



# CentriVap® Centrifugal Concentrators





# CentriVap® Required Accessories (not included)



Cold Trap

Vacuum  
Pump



Rotor



# CentriVap® Mobile Systems



## CentriVap Mobile Console

- Only true mobile system available
- Accessory vacuum pump fits inside cabinet
- Accommodates a variety of rotors

\*Exclusive feature

**LABCONCO**

# CentriVap® Mobile Systems



## CentriVap Mobile Console includes:


- Concentrator
- -55° C Cold Trap
- 12 - 13 mm Rotor
- Clear Canister with Solvent Insert
- New Vacuum Port
- Glass lid



# CentriVap® Centrifugal Concentrator Systems

Complete CentriVap Concentrator Systems are available for the following applications:

- Aqueous
- Solvent
- Acid
- Gel Dryer



\*All components necessary to begin operation are included except for the vacuum pump, which is sold separately.



# Gel Dryer Systems



CentriVap Gel Dryer Systems include the CentriVap Benchtop Concentrator, Cold Trap, Gel Dryer and connecting tubing.

The diaphragm pump shown above is sold separately.



# CentriVap® Centrifugal Concentrator Systems

## CentriVap Benchtop System includes:

- -55° Cold Trap
- 12 - 13 mm Rotor
- Clear Canister
- Chemical Insert\*—Moisture, Solvent or Acid



\* With Acid System, a glass insert will be included to protect Cold Trap





# CentriVap® Centrifugal Concentrator Systems


*Ultra-Low Cold Trap*



## Options & Accessories

- Ultra-Low Cold Trap (-85°C)
- Three-Way Valve allows connection of other laboratory equipment to CentriVap System
- Gel Dryer

*Gel Dryer*



*Three-way Valve*

# CentriVap® Cold Traps

## Cold Trap



- Epoxy-coated steel housing
- CFC-free refrigeration system with 1/4 hp compressor lowers cold trap to  $-55^{\circ}\text{C}$  in less than 15 minutes
- Stainless steel collector
- Ready indicator illuminates when collector coil reaches operating temperature of  $-40^{\circ}\text{C}$

*Lightbulb feature!*

# CentriVap® Cold Traps

## Ultra-Low Cold Trap

**Specs are same as the -50° Cold Trap with following exceptions:**

- Epoxy-coated steel housing and control panel
- CFC-free refrigeration system with one 1/4 hp and one 1/3 hp compressor lowers cold trap to -85° C in less than 30 minutes



# Acid Resistant CentriVap®



- Uses robust component materials to create corrosion resistant concentrator
- Resists exposure to TFA, HCl, Acetic Acids
- Includes Teflon-coated chamber, and selection of rotors, Hastelloy shaft, ceramic ball bearings and glass lid



# Refrigerated CentriVap® Centrifugal Concentrator

***Exclusive!***

**NEW!**



- For molecular biology applications with DNA, RNA and proteins
- Only concentrator of its kind!
- Heats to 100°C, cools to –4°C for heat sensitive samples
- Two-step time and temperature programming
- Pre-cooling option



# CentriVap® Rotors

Anodized aluminum rotors hold tubes at an angle to maximize surface area. Custom rotors can be ordered.

Acid resistant rotors available



Hexagonal Rotor



DNA Rotor

## CentriVap® Rotors

Rotors are interchangeable. The 12-13 mm Rotor is included with each CentriVap Aqueous, Acid, Solvent, Mobile and Gel Dryer System.



28 mm Rotor



12-17 mm Rotor



12-13 mm Rotor



# CentriVap® Rotors



**Microtiter Plate Rotor** hold four standard 96-well plates or two deep well plates



# CentriVap® Accessories



**2-Port Manifold**  
allows freeze  
drying of small  
volume samples in  
freeze dry  
glassware



**Glass Trap**  
Protects Cold  
Trap reservoir and  
lid from  
corrosives

**Clear  
Polypropylene  
Canister** holds  
chemical traps to  
provide additional  
protection to the  
vacuum pump



**Vacuum Controller**  
provides manual  
control and  
monitoring of  
vacuum levels



# Vacuum Pumps



**Diaphragm Vacuum Pumps** are Teflon-coated for corrosion resistance. Select pump with explosion-proof motor if using flammable solvents.



**Rotary Vane Direct Drive Vacuum Pumps** are two-stage, oil sealed and sliding-vane. A one liter bottle of pump oil and two inlet adapters are provided.



# Applications

| Type of Account   | Possible Application                         | Chemicals Used                 |
|-------------------|--|--------------------------------|
| Molecular Biology | Purification of Oligonucleotides             | Ammonium hydroxide<br>Aqueous  |
| Genetics          | Purification and concentrator of DNA and RNA | Methanol<br>Ethanol<br>Aqueous |
| Pharmaceutical    | Drug discovery<br>High throughput screening  | TFA<br>DMSO<br>DMF             |



# Norway - Chiron AS



## Kontakt oss

Tlf. 73 87 44 90 (Trondheim)

Tlf. 66 90 36 00 (Oslo)

Epost: [chiron@chiron.no](mailto:chiron@chiron.no)

Nettbutikk: [www.chiron.no](http://www.chiron.no)