



Product: cromatec

Art. No.: # 1700000

Commodity Code: 90189084

Highlights

- High throughput in routine H&E staining and cytology
- Continuous loading at any time to match individual workflows
- Individual timing station by station for utmost flexibility
- Enhanced Reagent Management System ensures best and consistent staining quality
- Fume extraction for your health and safety
- Easy operation via 7" color touch screen

Standard components

Basic instrument

| |
|--|
| 1 x Loading station |
| 1 x Unloading station |
| 2 x Slide racks for 30 slides each |
| 6 x Washing cuvettes (400 ml) |
| 20 x Staining cuvettes (400 ml) |
| 1 x Supply hose |
| 1 x Outlet hose |
| Fume extraction with 1 charcoal filter |
| Mains cable |
| Operation manual |
| Transport handles |
| Rack carrier including container |
| Adapter for most common racks |

Accessories

| | |
|----------------------------|-----------|
| Staining cuvettes (400 ml) | #41000017 |
| Washing cuvettes (400 ml) | #41000019 |
| Slide rack for 30 slides | #12000261 |
| Slide rack holder | #41000504 |
| Exchange carbon filter | #41000020 |

Technical Specification

General

| | |
|----------------------------------|--|
| Nominal voltage | 100 – 240 V AC +/- 10% |
| Nominal frequency | 50 / 60 Hz |
| Max. power draw | 250 VA |
| Max. heat emission | 250 J/s |
| IEC 1010 classification | Protective class 1 Pollution degree 2 Overvoltage installation category II |
| Fuses | 2 x T 2,5 A |
| Interfaces | USB, RJ45 |
| Operating temperature range | +10 to +35°C |
| Relative humidity | max. 80% non condensing |
| Temperature range during storage | +5°C bis +55°C |
| Humidity during storage | < 80% |

Dimensions and weight

| | |
|---------------------------------------|---------------------------|
| Dimensions (W x D x H) | 1560 mm x 640 mm x 620 mm |
| Weight unpacked (without accessories) | 110 kg |

Capacity

| | |
|--|-------------------------------------|
| Specimen slide throughput | Up to 1800 slides/h* |
| Loading capacity | Up to 20 slide racks simultaneously |
| Slide rack capacity | 30 slides |
| Total number of stations | 20 |
| Number of washing stations | Max. 6 |
| Reagent / Washing cuvettes volume | 400 ml |
| Water flow regulation for washing stations | 0,5 to 2,0 l/min |
| Load / unload stations | 5 each |
| Programs | 20 programs, up to 20 steps each |
| Incubation time setting | 0 sec. up to 59 min., 59 sec. |

* Depending on staining protocol