

Dear partner,  
we are pleased to inform you about the latest new features and developments for the Embedding center product line of SLEE medical.

slidetec  
WATER/HEAT

slidetec  
WATER

slidetec  
HEAT

The new histology water bath and slide warmer will add to the high quality of our products and increase your success in marketing and sales of SLEE medical products. The ergonomic design combined with precise temperature control allows a comfortable and productive processing of sectioned histological samples.

With Best Regards

SLEE medical GmbH

15.06.2015




## Maximize your process quality with highend instrumentation.

“Floaters” represent a potential cause of diagnostic error and may arise during tissue processing and slide preparation.

Contaminant tissue in microscopic slides is a potential source of diagnostic error. When such extraneous tissue is neoplastic, it may result in a false-positive diagnosis. In addition, misplaced fragments of malignant tissue originating from the patient’s specimen may result in upstaging.

Layfield et al., Am J Clin Pathol 2011;136:767-772



**The histology water bath slidetec<sub>WATER</sub> with integrated illumination (LED (white/blue)) allows the selection of the best sections and supports prevention of cross-contamination.**




## slidetec

### WATER/HEAT

The SLEE floating bath and slide warmer combination slidetec WATER/HEAT is designed for stretching of paraffin sections and subsequently drying and baking these on microslides. It combines functional design with safety and ergonomic features that improve process reliability and convenient use. With slidetec HEAT and slidetec WATER also the single units are available.

- **EASY WORK** Ergonomic oriented and spacious drying area for slides
- **EASY VIEW** Integrated water bath illumination
- **EASY CLEAN** Separable glass water bath
- **EASY SAFE** Overheat protection
- **EASY PRECISION** Microprocessor thermostat control
- **EASY EFFICIENCY** Programmable Working hours



### Specifications

|   |                    |
|---|--------------------|
| <b>Capacity water bath</b>                | 2 liters           |
| <b>Capacity slide warmer</b>              | 30 slides          |
| <b>Dimension (width x depth x height)</b> | 560 x 370 x 230 mm |
| <b>Weight</b>                             | 23 kg              |
| <b>Temperature range</b>                  | RT to 60 °C        |

|  |                           |
|--|---------------------------|
| <b>slidetec WATER/HEAT 230 V, 50/60 Hz</b> | <b>Art. No. #13000500</b> |
| <b>slidetec WATER/HEAT 115 V, 50/60 Hz</b> | <b>Art. No. #13000501</b> |



## slidetec WATER

### Specifications

|   |                    |
|---|--------------------|
| <b>Capacity water bath</b>                | 2 liters           |
| <b>Dimension (width x depth x height)</b> | 350 x 370 x 230 mm |
| <b>Weight</b>                             | 10 kg              |
| <b>Temperature range</b>                  | RT to 60 °C        |

|                                       |                           |
|---------------------------------------|---------------------------|
| <b>slidetec WATER 230 V, 50/60 Hz</b> | <b>Art. No. #13000520</b> |
| <b>slidetec WATER 115 V, 50/60 Hz</b> | <b>Art. No. #13000521</b> |



## slidetec HEAT

### Specifications

|   |                    |
|---|--------------------|
| <b>Capacity slide warmer</b>              | 30 slides          |
| <b>Dimension (width x depth x height)</b> | 350 x 370 x 230 mm |
| <b>Weight</b>                             | 15 kg              |
| <b>Temperature range</b>                  | RT to 60 °C        |

|                                      |                           |
|--------------------------------------|---------------------------|
| <b>slidetec HEAT 230 V, 50/60 Hz</b> | <b>Art. No. #13000510</b> |
| <b>slidetec HEAT 115 V, 50/60 Hz</b> | <b>Art. No. #13000511</b> |

## **INSTRUMENTS**

slidetec HEAT/WATER 230 V, 50/60 Hz

Art.No.: #13000500

slidetec HEAT/WATER 115 V, 50/60 Hz

Art.No.: #13000501

slidetec HEAT 230 V, 50/60 Hz

Art.No.: #13000510

slidetec HEAT 115 V, 50/60 Hz

Art.No.: #13000511

slidetec WATER 230 V, 50/60 Hz

Art.No.: #13000520

slidetec WATER 115 V, 50/60 Hz

Art.No.: #13000521

SLEE medical GmbH

Lise-Meitner-Strasse 11  
55129 Mainz · Germany

Tel. +49 6131 95871-0  
Fax +49 6131 95871-722

[mail@slee.de](mailto:mail@slee.de)  
[www.slee.de](http://www.slee.de)