

Multiple functions for easy and reliable pipetting!



Transferpette® α

F I R S T C L A S S · B R A N D

The Transferpette® α pipettes have all important features essential for demanding Life Science applications.

■ **Ergonomic finger rest**

Perfectly adapted to the shape of the hand

■ **Clearly visible volume display**

4 digits for precise volume setting. Always visible, even while pipetting – for right- and left-handed users.

■ **Autoclavable without disassembling**

The Transferpette® α is completely autoclavable!

■ **Full range from 0.1 μ l up to 10 ml**

The entire range can be covered with only 5 models.



Transferpette® α – light, high-precision, robust and dependable

Particular advantages include the convenient adjustment and operation, 4-digit volume setting and full autoclavability.

The piston operates smoothly, for effortless pipetting with a light touch. All of these features make the Transferpette® α a convenient, very reliable pipette in everyday lab routine.



- Corrosion-resistant piston
- Innovative plastic materials
- Centrally located pipetting button and separate ejection key
- CE-IVD-certified
- Transferpette® α 0.1 - 2.5 μ l – highest precision for molecular biology, especially for pipetting enzymes.

Technical Data

Items supplied:

Transferpette® α , adjustable volume, DE-M marking, performance certificate, silicone oil.


Capacity μ l	Description	Type of tip μ l	Nominal volume μ l	A* $\leq \pm$ %	CV* \leq %	Subdivision μ l	Cat. No.
0.1 - 2.5	D-2.5 ③	20 / 50 / nano-cap™	2.5	1.4	0.7	0.002	7049 69
0.5 - 10	D-10	20 / 50 / nano-cap™	10	1	0.5	0.01	7049 70
2 - 20	D-20	200	20	0.8	0.4	0.02	7049 72
5 - 50	D-50	200	50	0.8	0.3	0.05	7049 73
10 - 100	D-100	200	100	0.6	0.2	0.1	7049 74
20 - 200	D-200	200	200	0.6	0.2	0.2	7049 78
100 - 1000	D-1000	1000	1000	0.6	0.2	1	7049 80
500 - 5000	D-5000 ①	5000	5000	0.6	0.2	5	7049 82
1000 - 10000	D-10000 ②	10000	10000	0.6	0.2	10	7049 84

Note:

① D-5000 for use only with 5 ml tip from BRAND and 5 ml tip from Thermo Fisher Scientific Inc. (Finnpipette®).

② D-10000 for use with 10 ml tip from BRAND and appropriate tips from Gilson® and Eppendorf®.

③ D-2.5: when using BRAND nano-cap™ tips, 50 μ l tips or tips from other manufacturers, calibration and adjustment are possibly necessary.



Shelf/rack mount for single Transferpette® α pipette.

Pack of 1.

Cat. No. 7048 10



Filter

for Transferpette® α 0.5-5 ml pipettes. Pack of 25.

Cat. No. 7046 52

for Transferpette® α 1-10 ml pipettes. Pack of 25.

Cat. No. 7046 53

Transferpette® and BRAND® are trademarks of BRAND GMBH + CO KG, Germany.

Other reproduced brands are the property of the respective owner.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG · P.O. Box 1155 · 97861 Wertheim · Germany

Tel.: +49 9342 808-0 · Fax: +49 9342 808-98000 · E-Mail: info@brand.de · Internet: www.brand.de

