

Multifarious Possibilities

MADE IN
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

Reax top

Powerful mixing – space saving – broadly usable

Support plate made of non-slip rubber enables safe vortexing for sample vessels with a Ø 20 mm

The support plate can be easily removed and replaced when dirty or broken

The change for the large test tube tray or attachment for 10 reaction vessels enables use for many applications

Rubber feet ensure a secure hold at high speeds and demanding mixing tasks

Thread for attaching the vessel holder. Test tubes can be shaken safely even in continuous operation

The slice switch allows you to choose between short-term and continuous operation

The speed controller enables comfortable operation



Accessoires



Large test tube tray for round bottom vessel up to 50 ml



Attachment for 10 test tubes



Test tube holding device



Test tube stand