

פיפטורים איכותיים

מסדרת Transferpette α החדשה,
תוצרת BRAND, גרמניה!



כל הפיפטור עמיד לעיקור באוטוקלאב, ללא צורך בפירוק.

ידיית אחיזה ארגונומית, להתאמה מושלמת לכף היד והפעלה במינימום כוח.

מפחית את הסיכון לדלקות בשורש כף היד (RSI).

הפעלה קלה ביד אחת, לימניים ושמאליים כאחד.

תצוגת הנפח ניכרת גם בעת השימוש.

כיוון נפח מדויק ביותר, מוצג ב-4 ספרות עשרוניות.

בוכנה עמידה במיוחד לחומרים קורזיביים.

עמיד בפני חיסוי ב-UV.

כל פיפטור מסומן בצבע בהתאם לטיפ הנדרש.

טווח נפחים רחב: 0.1ul-10ml


Capacity μ l	Description	Type of tip μ l	Nominal volume μ l	A* $\leq \pm$ %	CV* \leq %	Subdivision μ l	Cat. No.
0.1 - 2.5	D-2.5 ①	20 /50 / nano-cap™	2.5	1.4	0.7	0.002	7049 69
0.5 - 10	D-10	20 /50 / nano-cap™	10	1	0.5	0.01	7049 70
2 - 20	D-20	200	20	0.8	0.4	0.02	7049 72
5 - 50	D-50	200	50	0.8	0.3	0.05	7049 73
10 - 100	D-100	200	100	0.6	0.2	0.1	7049 74
20 - 200	D-200	200	200	0.6	0.2	0.2	7049 78
100 - 1000	D-1000	1000	1000	0.6	0.2	1	7049 80
500 - 5000	D-5000 ②	5000	5000	0.6	0.2	5	7049 82
1000 - 10000	D-10000 ③	10000	10000	0.6	0.2	10	7049 84

Note:

① D-5000 for use only with 5 ml tip from BRAND and 5 ml tip from Thermo Fisher Scientific Inc. (Finnpipette®).

② D-10000 for use with 10 ml tip from BRAND and appropriate tips from Gilson® and Eppendorf®.

③ D-2.5: when using BRAND nano-cap™ tips, 50 μ l tips or tips from other manufacturers, calibration and adjustment are possibly necessary.



לבוטל ציוד מדעי בע"מ

לפרטים נוספים והדגמות: